

Portable Rice Husker



SCIENCE OF SENSING

Model TR130



Tomorrow's Testing, Today !
1-800-438-5388

Portable Rice Husker

Easily Husk Small Rice Samples

This simple-to-use instrument husks small paddy rice samples to provide brown rice. This process is used to generate samples for accurate moisture measurements as well as for many unique R&D applications.

Simple Operation - World-Class Results

To husk a sample, simply pour paddy rice into the TR130. Close the lid and turn the crank. While husking the rice, the brown rice is automatically moved to a separate case where it can be removed for the moisture test. The husk can be easily discarded.

The operation is very clean as the test is completely enclosed. Each sample is approximately 8g. Since manual sample preparation is eliminated, using the TR130 improves the efficiency of a rice testing laboratory or grain elevator as it makes it easy to take a variety of moisture measurements. In addition, test results are more stable because of the improved brown rice consistency.

Other Kett Rice Testing and Processing Instruments

The TR130 is the latest design from Kett, considered the world leader in rice quality measurement. We manufacture a wide array of instruments and processing aids. Some of our more popular products are shown below.

Specifications

Model	TR130
Principle	Grinding, Husking
Sample Volume	Approximately 8g
Product	Paddy Rice (to Brown Rice)
Moisture Range	12-18%
Husking Time	Approximately 20 seconds
Weight (kg) (Net/Shipment)	0.2/1
Dimensions (mm)	112(W)x79(H)x70(D)
Warranty	One Year Total Care Maintenance Agreements available



• PM450 Portable
Moisture Meter



• PM650 Advanced
Moisture Meter



• Riceter portable
Moisture Meter



• AN820 Instant Rice
Composition Analyzer



• C600 Rice
Whiteness Meter

Built to Last - Lab Quality Results - Production Line Durability

Kett stands behind the TR130 systems with a full **One-Year Total Care** manufacturer's warranty on parts, labor and specification compliance. Kett's dependability and reliability allow you to productively focus on improving your product's quality and not waste time calibrating and verifying our gauges. Kett is recognized as a worldwide leader in grain measurement instrumentation. With over 65 years of design and manufacturing excellence, the TR130, like all Kett instruments, is the standard by which others are measured.



• TX200 Rice Analyzer



• RN300 Rice
Quality Analyzer



• RN600 Automated
Quality Analyzer



• TR250 Automatic
Rice Husker



• Pearlest Rice Polisher

Tomorrow's Testing, Today!

In addition to our unique TR130, Kett manufactures almost 200 different instruments. Available for portable handheld, laboratory desktop or online process applications, if you have a measurement requirement, we offer solutions...today! If you need to measure organic constituents - - accurately, quickly, and simply, wish to place an order, or require additional guidance on model selection, please call Kett toll-free!

1-800-GET-KETT



PMB 504, 17853 Santiago Boulevard,
Suite 107, Villa Park, CA 92861
714-974-8837 • 714-974-8877 (FAX)

1-800-438-5388

www.Kett.com

Sales, Support, Service

ISO 9001

